



Packaging Solutions As Uniquely Designed As The Devices You Develop



## HEADER BAGS

Medical Grade Header Bags Ensure Aseptic Presentation for Medical Device, Surgical Kits and Pharmaceutical Applications Requiring ETO or Radiation Sterilization.



Techni-Vent Header Bags are made of DuPont™ Tyvek® for strength and bacterial barrier, ensuring clean visual transfer of coating which confirms pouch integrity and peelability with no fiber tear or particulate. They are made of puncture and abrasion resistant film and laminates under strict ISO 13485 production standards. Our header bags can be manufactured to accommodate irregular and/or high-profile items.



[www.technipaq.com](http://www.technipaq.com)

975 Lutter Drive  
Crystal Lake, IL 60014

Phone: 815.477.1800  
Fax: 815.477.0777  
[info@technipaq.com](mailto:info@technipaq.com)

LAMINATING

COATING

PRINTING

SLITTING

CONVERTING



MEDICAL  
DEVICE



DIAGNOSTIC



LIFE  
SCIENCE



PHARMACEUTICAL

## TO MEET YOUR SPECIFIC REQUIREMENTS

We offer a wide range of films, laminates, and foil structures in a variety of sizes, as well as DuPont™ Tyvek® Medical & Pharmaceutical Packaging materials Tyvek® 1073B, Tyvek® 1059B, Tyvek® 2FS™ and Tyvek® 40L.

### CAPABILITIES

- Stock sizes available for immediate shipment – with Tyvek® vents
- Small or large quantity runs
- Unprinted or printed pouches in up to eight colors
- Custom pouches available up to 48"x99"

### MATERIALS

- DuPont™ Tyvek® Nylon, polyethylene, or high puncture resistance options
- Specialty films
- Peelable laminates
- Foil laminates
- Clear high barrier

### STERILIZATION

- Radiation
- Autoclave
- ETO

### DUPONT™ TYVEK® MATERIALS FEATURES AND BENEFITS

- Outstanding resistance to microbial penetration
- Superior tear strength and puncture resistance
- Compatible with a broad range of sterilization methods

### STOCKED SIZES

Featuring 1073B DuPont™ Tyvek® with 4mil LLDPE

8" x 12" w/3.75"

10" x 15" w/3.75"

12" x 16" w/3.75"

13" x 15" w/3.75"

16" x 22" w/3.75"

19" x 30" w/3.75"

**We do  
what others  
don't.**

**We customize  
where others  
won't.**



**DESIGN YOUR OWN  
POUCH IN MINUTES**  
[technipaq.com/design](http://technipaq.com/design)

«DUPONT»  
**Tyvek**®



[www.technipaq.com](http://www.technipaq.com)

975 Lutter Drive  
Crystal Lake, IL 60014

Phone: 815.477.1800  
Fax: 815.477.0777  
[info@technipaq.com](mailto:info@technipaq.com)